BPL_CHO_Perfusion_cspr_openloop script with FMPy

The key library FMPy is installed.

After the installation a small application BPL_CHO_Fedbatch is loaded and run. You can continue with this example if you like.

!lsb_release -a # Actual VM Ubuntu version used by Google No LSB modules are available. Distributor ID: Ubuntu Description: Ubuntu 22.04.3 LTS 22.04 Release: Codename: jammy %env PYTH0NPATH= → env: PYTH0NPATH= !wget https://repo.anaconda.com/miniconda/Miniconda3-py312_24.3.0-0-Linux-x86_64. !chmod +x Miniconda3-py312_24.3.0-0-Linux-x86_64.sh !bash ./Miniconda3-py312 24.3.0-0-Linux-x86 64.sh -b -f -p /usr/local import sys sys.path.append('/usr/local/lib/python3.12/site-packages/') → --2024-05-23 09:52:27-- https://repo.anaconda.com/miniconda/Miniconda3-py312 Resolving repolanaconda.com (repolanaconda.com)... 104.16.191.158, 104.16.32.2 Connecting to repo.anaconda.com (repo.anaconda.com)|104.16.191.158|:443... com HTTP request sent, awaiting response... 200 OK Length: 143351488 (137M) [application/octet-stream] Saving to: 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh' Miniconda3-py312_24 100%[===========] 136.71M 31.7MB/s in 4.5s 2024-05-23 09:52:32 (30.1 MB/s) - 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh' PREFIX=/usr/local Unpacking payload ... Installing base environment... Preparing transaction: ...working... done Executing transaction: ...working... done installation finished. !conda update -n base -c defaults conda --yes → Channels: - defaults Platform: linux-64

Collecting package metadata (repodata.json): done

Solving environment: done

Package Plan

environment location: /usr/local

added / updated specs:

- conda

The following packages will be downloaded:

| package | build | |
|----------------------------------------------------|--------------------------------------------------|---------------------------|
| conda-24.5.0 frozendict-2.4.2 openssl-3.0.13 | py312h06a4308_0 py312h06a4308_0 h7f8727e_2 | 1.2 MB 36 KB 5.2 MB |
| | Total: | 6.5 MB |

The following NEW packages will be INSTALLED:

frozendict pkgs/main/linux-64::frozendict-2.4.2-py312h06a4308_0

The following packages will be UPDATED:

conda 24.3.0-py312h06a4308_0 --> 24.5.0-py312h0 openssl 3.0.13-h7f8727e_0 --> 3.0.13-h7f872

Downloading and Extracting Packages:

openssl-3.0.13 | 5.2 MB | : 0% 0/1 [00:00<?, ?it/s] conda-24.5.0 | 1.2 MB | : 0% 0/1 [00:00<?, ?it/s]

openssl-3.0.13 | 5.2 MB | : 0% 0.002992912657148922/1 [00:00<00:48

frozendict-2.4.2 | 36 KB | : 44% 0.43853215920344746/1 [00:00<00:00, conda-24.5.0 | 1.2 MB | : 1% 0.01293349794914382/1 [00:00<00:11,

frozendict-2.4.2 | 36 KB | : 100% 1.0/1 [00:00<00:00, 2.89it/s] conda-24.5.0 | 1.2 MB | : 100% 1.0/1 [00:00<00:00, 1.87it/s]

Preparing transaction: done Verifying transaction: done Executing transaction: done

!conda --version
!python --version

conda 24.5.0 Python 3.12.2

!conda install -c conda-forge fmpy --yes # Install the key package



Preparing transaction: done Verifying transaction: done Executing transaction: done

```
#!conda install matplotlib --yes

#!conda install scipy --yes

#!conda install openpyxl --yes

#!conda install xlrd --yes
```

BPL_CHO_Perfusion_cspr_openloop setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:

- FMU BPL_CHO_Perfusion_cspr_openloop_linux_om_me.fmu
- Setup-file BPL_CHO_Perfusion_cspr_openloop_fmpy_explore

```
# Filter out DepracationWarnings for 'np.float as alias' is needed — wish I could import warnings warnings.filterwarnings("ignore")

%%bash git clone https://github.com/janpeter19/BPL_CHO_Perfusion

→ Cloning into 'BPL_CHO_Perfusion'...

%cd BPL_CHO_Perfusion

→ /content/BPL_CHO_Perfusion
```

BPL_CHO_Perfusion_cspr_openloop - demo

Author: Jan Peter Axelson

```
run -i BPL_CHO_perfusion_cspr_openloop_explore_fmpy.py

Linux - run FMU pre-comiled OpenModelica 1.23.0-dev

Model for bioreactor has been setup. Key commands:
- par() - change of parameters and initial values
- init() - change initial values only
- simu() - simulate and plot
- newplot() - make a new plot
- show() - show plot from previous simulation
- disp() - display parameters and initial values from the last simulation
- describe() - describe culture, broth, parameters, variables with values/un
```

Note that both disp() and describe() takes values from the last simulation and the command process_diagram() brings up the main configuration

Brief information about a command by help(), eg help(simu)
Key system information is listed with the command system_info()


```
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
```

1 About the process model

We can get information about the process and liquid phase by the command describe(). Here is no gas-phase included. This command can also be used to bring up information about a specific variable or parameter. However, you should use describe() after a simulation to get the valued used during the simulation.

```
process_diagram()
```

No processDiagram.png file in the FMU, but try the file on disk.



describe('culture'); print(); #describe('liquidphase')

Reactor culture CHO-MAb - cell line HB-58 American Culture Collection ATCC


The molecular weight of the recombinant protein (MAb) is somewhat arbitrarly chosen and the value not used in the simulations.

describe('MSL')

→ MSL: 3.2.3 – used components: RealInput, RealOutput, CombiTimeTable, Types

newplot('Perfusion cultivation', plotType='Extended')

```
par(samplePeriod=1); par(F1=0.0020)  # General parameters
par(mu_ref=0.04); simu()  # First simulation (solid)
par(mu_ref=1); simu()  # Second simulation (dashed)
```




model_get('CSPR')

→ 105.41563602162672

We see here the shorter settling time for cell conc Xv when the slower exponential increase of perfusion rate is used, compared to an abrupt change.

newplot('Perfusion cultivation', plotType='Extended')


```
par(samplePeriod=1); par(F1=0.0080)  # General parameters
par(mu_ref=0.04); simu()  # First simulation (solid)
par(mu_ref=1); simu()  # Second simulation (dashed)
```



We see that a somewhat smaller change in perfusion rate make the differnce in settling time of Xv for the two strategies to almost disappear.

```
newplot('Perfusion cultivation', plotType='Extended')
simu(600)
simu(400,'cont')
```





describe('feedtank.W')

```
→ Variable not logged
describe('feedtank.V')
Feed volume : 78.577 [ L ]
    ---- I MESCLINE ( LEGALATIVIM )
describe('bioreactor.broth_decay.k_decay')
Rate of decay of viable cells : 0.0 [ 1E9 cells/(L*h) ]
                except AttributeError:
describe('k_lysis')
→ Specific rate of lysis of dead cells : 0.0 [ 1/h ]
disp('culture')
\rightarrow qG_max1 : 0.297
    qG_max2 : 0.038
    qGn_max1 : 0.124
    qGn_max2 : 0.022
    mu_d_max : 0.13
describe('mu')
Specific cell growth rate variable: 0.026 [ 1/h ]
describe('parts')
['bioreactor', 'bioreactor.broth_decay', 'bioreactor.culture', 'CSPR', 'cspr_o
```